







LSFSF98UK

UNDER BENCH LABORATORY FREEZER
SOLID DOOR 98L WHITE



- 
-18°C to -20°C Operating Temperature Range
 Conforms to MHRA guidelines.
- 
Lock Fitted as Standard
 Maintains security and ensures the unit can only be opened by authorised personnel.
- 
Unique Anti Microbial Handle
 Helps prevent the growth of bacteria and minimise the potential for cross contamination
- 
3 Year Extended Warranty
 When serviced by an approved service technician, the warranty is extended to 3 years.



The Pharmacy
Guild of Australia



Specifications

Capacity	98L Gross / 84L Net	Spark Free	Yes
Dimensions External	595w x 660d x 850h	Internal Fan	No
Dimensions Internal	420w x 430d x 640h	Temperature range	-18°C – +20°C
Lock	Yes, 2 Keys	Controller	Digital, External
Door	Solid	Interior	HIPS
Door Style	Right Hand Hinge (Reversible)	Housing	Powder coated sheet
Weight	37.25 kg	Voltage 50/60Hz	230V
Drawers fit/max	3	Internal light	No
Drawers Material	Plastic	Door Open Alarm	Yes
Defrost	Manual	Access Port	No
Refrigerant	R600A	Noise Level (db)	-
Adjustable Feet (No & Position)	4 adjustable feet	Anti Microbial Handle	Yes
Country of Origin	United Kingdom	Colour	White / Green

ROLLEX MEDICAL NZ LTD

PO Box 204159, Highbrook, Auckland 2161, New Zealand
 Phone: +64 9 274 0707 Fax: +64 9 274 0709
 Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

ROLLEX MEDICAL (AUST) PTY LTD

PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia
 Phone: 1300 880 441 | Fax: 1330 880 451
 Email: sales@rollexmedical.com | www.rollexmedical.com.au